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
Quantitative Imaging of Trace Elements in Mouse Spinal Cord using fsLA-ICP-MS

Femtosecond (fs) LA-ICP-MS is a powerful technique for the trace element quantification in hydrous samples due to the under atmospheric pressure measurement. Moreover, isolated sampling and ionization processes, which eliminate the influence of coexisting component, are capable of high sensitivity and accuracy.

Principle of fsLA-ICP-MS


*) Courtesy of Prof. Hirata, University of Tokyo.

Advantages of fsLA-ICP-MS

	fsLA-ICP-MS	NanoSIMS	TOF-SIMS
Sensitivity	ppb – sub ppm	ppm	ppm
Spatial Resolution	10 μ m	50 nm	300 nm
Atmosphere	Atmospheric pressure, He (Available for Hydrous Samples)	Ultra high vacuum (< 10^{-7} Pa)	Ultra high vacuum (10^{-6} – 10^{-7} Pa)
Sampling Area	~mm – 2 cm	10 – 50 μ m	10 μ m – 5 cm
Sampling Depth	~50 μ m	< 10 nm	< 3 nm

- Trace elemental imaging with **High Sensitivity** in **Biological Tissue**
- Quantitative analysis** with isolated sampling and detection region (Less influence of coexisting components)

Homogenized Biological Standard Preparation & Calibration



Calibration for Spinal Cord

- 1st : Fe/Zn concentrations of ③ determined by ICP-MS, Fe/Zn/C Intensities of ④ measured by fsLA-ICP-MS.
- 2nd : Matrix corrections between brain and spinal cord are provided on the basis of carbon concentrations.

- Standard with Targeted Concentration can be prepared.

Quantitative Imaging of Fe/Zn in Lumbar Cord

- Fe/Zn are more highly concentrated in **Gray matter** (10–15 ppm) than **White matter** (3–5 ppm).

fsLA-ICP-MS & reliable standard preparation make **Quantitative Element Imaging** possible with **High Sensitivity** and **Accuracy**.


Candidate method for **Drug Efficacy Assessment**, **Pathological Mechanisms** analyses, etc.

SEC-MALS Analysis under GMP — Molecular Weight Determination of Aggregates —

Evaluation of aggregates is indispensable for the quality test of antibody drugs. Absolute molecular weight distribution, mean absolute molecular weight, and radius of gyration of macromolecules are determined by SEC-MALS without constructing a calibration curve. In this report, aggregates were prepared by heating the human IgG antibody, and the molecular weight of the aggregates was determined by SEC-MALS.

Principle of SEC-MALS

[MALS]


Multi Angle Light Scattering (MALS) detector:
Scattered light is detected by 18 photodiodes.

Refractive Index (RI) detector for measuring concentration

$$\frac{K^*C}{R(\theta)} = \frac{16\pi^2 n_0^2}{3M\lambda_0^2} < r_g^2 > \sin^2\left(\frac{\theta}{2}\right) + \frac{1}{Mw}$$


Mean absolute molecular weight (Mw) and radius of gyration (r_g) are calculated from concentration (measured by RI) and light scattering intensity (measured by MALS).


Zimm plot

We have introduced a SEC-MALS system administrated under GMP as a Japanese CRO for the first time. The acquired data are available for the New Drug Application.


Analysis of BSA (For MALS calibration)


Analysis of Antibody

Analysis of Molecular Weight Markers

Analysis of Antibody Aggregates

GC/MS Analysis under GMP – Highly Sensitive Analysis for Impurities in Nucleic Acids

In the quality test of nucleic acid medicines, plenty of samples can not be used for analysis, unlike the case of small molecule drugs. Especially in the analysis for impurities, highly sensitive analysis is needed to perform the test using limited amount of valuable samples. In this report, the residual solvents contained in a commercially available nucleic acid sample were measured by highly sensitive headspace-GC/MS.

1. Quality Test Items for Nucleic Acid Medicines


Requirement of guideline*	Available analysis at Toray Research Center
Specification: Test procedures	
(1) Appearance	Observation
(2) Identification	UHPLC, UV spectrophotometer
(3) Purity	
1. Oligonucleotide related substance	UHPLC, Capillary electrophoresis
2. Organic small molecular impurities	LC-MS/MS, GC-MS, MS, NMR, UHPLC
3. Residual solvents	GC-MS, GC
4. Inorganic impurities	ICP-OES, ICP-MS
(5) Assay	UHPLC, UV spectrophotometer
(6) Biological activity	Cell-based assay, ELISA

*Quoted from "Matters to be considered for assurance and evaluation of the quality of nucleic acid medicines" issued by MHLW, Japan on September 27, 2018


SOP: Operation and maintenance of headspace-GC/MS (SOP No. NME18500)

2. Measurement of Residual Solvents in the Nucleic Acid by Headspace-GC/MS

<Figure of a Headspace Device>


Improvement of detection capability for dichloromethane (a typical solvent used for nucleic acid synthesis)

Chromatograms of the dichloromethane solution

Residual dichloromethane in artificially synthesized nucleic acid was measured by the trap mode.

Sample: Artificially synthesized nucleic acid (C-Myc antisense DNA, sulfurized)
Sequence: CACGTTGAGGGGCAT
Sample concentration: 1 mg/mL

Chromatogram of the artificially synthesized nucleic acid by the trap mode

A trace amount of dichloromethane (4ppm for the synthesized nucleic acid) was detected.

- S/N ratio of dichloromethane increased by 30 times in measurement by trap mode headspace-GC/MS.
- Quantitation of dichloromethane was performed using very small amount (1 mg) of the nucleic acid.

Analysis of amino acids, metabolic and buffer constituents in a medium

A medium contains low molecular weight nutrients, balanced saline and buffer constituents, which can be quantified by means of capillary electrophoresis (CE). By applying diverse separation conditions, comprehensive analysis of ionic species, such as HEPES commonly used for a buffer and amino acids required to control for cell culture, in a medium can be achieved.

I. Ionic species in a medium

Examples of evaluable ionic species

The colors of letters (red, green, blue, purple) are corresponding to A, B, C, D in the section II.

Amino acids

Glutamine (Gln), Glutamate (Glu)
 Pyroglutamate, Glycine (Gly)
 Histidine (His), Methionine (Met)
 Arginine (Arg), Lysine (Lys)
 etc.

Metabolic constituents

Pyruvate, Lactate
 Citrate, Malate
 Succinate
 Choline, Carnitine
 etc.

Buffer constituents

Acetate, PO_4^{3-}
 Citrate, Succinate
 HEPES
 etc.

Minerals


Cl^- , PO_4^{3-} , SO_4^{2-}
 NO_3^- , Na^+ , Ca^{2+}
 K^+ , Mg^{2+} , NH_4^+
 etc.

Media


II. CE analyses of amino acids, metabolic and buffer constituents in media

(A) ✓ Inorganic anions, metabolic and buffer constituents
 ✓ Acidic amino acids like pyroglutamate


(B) ✓ Amino acids, HEPES*
 ✓ Pyroglutamate, Glu, Gln

*HEPES : 2-[4-(2-Hydroxyethyl)-1-piperazinyl]ethansulfonate


Amino acids for nutrients, as well as inorganic and organic anions as metabolic and buffer constituents, can be evaluated.

(C) ✓ Inorganic cations, metabolic constituents
 ✓ Basic amino acids like Arg

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Inorganic cations, amines and basic amino acids as metabolic constituents can be evaluated.

Fig. Electropherogram of ionic species in media (A), (B) anionic species (C), (D) cationic species

CE based on optimized conditions enables the evaluation of a lot of ionic species with the LOQ** 0.5 mg/mL.

**LOQ : limit of quantification